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### Giriş

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## References:

1. S. Pors Nielsen, "The Fallacy of BMD: A Critical Review of the Diagnostic Use of Dual X-Ray Absorptiometry," *Clinical Rheumatology* 2000 19:174-183
2. Ott S. M. et al., "Evaluation of Vertebral Volumetric vs. Areal Bone Mineral Density During Growth," *Bone* 1997 June 20(6):553-6
3. Pors Nielsen et al., "Bone Densitometry – Two or Three Dimensions?" In: *Current Research in Osteoporosis and Bone Mineral Measurement V*, Ring E. F. J., Elvins D. M. and Bhalla A.K. (eds) pp. 34-35. London: British Institute of Radiology.
4. Formica C. A., "Standardization of BMD Measurements," Editorial, *Osteoporos International* 1998; 8:1-3.
5. Blake G. M. et al., "An Unexpected Change in DXA Calibration Not Detected by Routine Quality Control Checks," *Osteoporos International* 1999; 9:115-120.
6. Kanis et al., "Guidelines for Diagnosis and Management of Osteoporosis," *Osteoporos International* 1997; 7:390-406.
7. Sievanen H., and the Bone Research Group at the UKK Institute in Tampere, Finland, "QUSDerived Speed of Sound and Cortical Bone Structure," (abstract) presented at the ASBMR 21st Annual Meeting in St. Louis, MI, USA, September 1999.

8. Bouxsein, M.L. et al., "Prediction of Femoral Failure Load from Femoral BMD and Ultrasonic Velocity at the Femur, Radius and Phalanx," (abstract), presented at the ASBMR 21 st Annual Meeting in St. Louis, MI, USA, September 1999.
9. Weiss M., Ben Shlomo A., Hagag P., Rapoport M., "Reference Database for Bone Speed of Sound Measurement by a Novel Quantitative Multi-site Ultrasound Device," *Osteoporosis International* 2000 11:688-696.
10. Knapp et al., "Multiple Site Ultrasound Measurements Predict Vertebral Fractures in Postmenopausal Women," (abstract), presented at the ASBMR 21st Annual Meeting in St.Louis, MI, USA, September 1999.
11. Knapp et. al., "Ultrasound Measurements at the Radius Predict Wrist Fractures in Postmenopausal Women," (abstract), presented at the Annual Meeting of the Bone and Tooth Society in Bristol, UK, June 1999.
12. Sunlight Omnisense™ FDA PMA (Pre-Market Application) Approval, January 20, 2000.
13. Weiss M. et al., "The Importance of Precision – New Hopes for Monitoring Osteoporosis Treatment by QUS," (abstract), presented at the ISCD meeting, Brazil 2000.
14. Weiss M., Ben Shlomo A., Hagag P., Ish-Shalom S., "Discrimination of Proximal Hip Fracture Measurement at the Radius," *Osteoporosis International* 2000 11:411-416.
15. Barkmann R. et al., "A New Method for Quantitative Ultrasound Measurements at Multiple Skeletal Sites – First Results of Precision and Fracture Discrimination," *Journal of Clinical Densitometry* 2000 spring (1) 1-7.
16. Knapp, K. et al., "Quantitative Ultrasound Measurements Detect Skeletal Changes in Cortical Bone Following HRT Use," (abstract), presented at the 11th International Workshop on Calcified Tissues, Eilat, Israel, February 1999.
17. Weiss M. et al., "Effect of Estrogen Replacement Therapy on Speed of Sound at Multiple Skeletal Sites," *Maturitas* 2000 (35) 237-243.
18. Ben Shlomo A., Weiss, M. et al., "Thyroid Dysfunctional State Detected by QUS Measurement at Multiple Skeletal Sites," (abstract), presented at the ASBMR 2nd Joint Meeting, California, USA, December, 1998. Also In Press.
19. M. Weiss et al., "Early Affect of Alendronate or Raloxifene Treatment in Osteoporotic Women Monitored by Multi-Site QUS," (Abstract) presented at the ASBMR 22 nd Annual Meeting in Toronto, Canada, September 2000.
20. Njeh C.F. et al., "An In Vitro Investigation of the Dependence on Sample Thickness of the Speed of Sound along the Specimen," *Medical Engineering & Physics*, 1999 September, 21: 651-659.
21. Hans, D. et al., "Does Combining the Results from Multiple Bone Sites Measured by a New Quantitative Ultrasound Device Improve Discrimination of Hip Fracture?" *Journal of Bone Mineral Research*, 1999, Vol. 14, No. 4 644-6511
22. Knapp, K. et al., "Preliminary Results Of The Sunlight Omnisense™ Bone Sonometer: In-Vivo and In-Vitro Precision and Correlation With DXA," (abstract), The ASBMR-IBMS 2 nd Joint Meeting, California, USA, December 1998.