

Measurement Report

7/2/99

Patient Information

First Name Helen
 Family Name Jackson
 Age 54y9m
 Height 5'5"
 Weight 132 lbs.
 BMI 22

ID..... 38756348
 Gender..... Female
 Referring Physician Dr. James Smith
 Reporting Physician Dr. Bernard Costner
 Operator Lynn Moore


Measurement Results

Site..... **DISTAL RADIUS**

SOS..... **3982[m/sec]**
 Speed of Sound

T-Score..... **-1.9**
 Number of standard deviations from the mean of healthy young adults

Z-Score..... **-1.1**
 Number of standard deviations from the mean of age matched population



Comments and Findings

Comments/Anamnesis:

Findings:

Normal

Osteopenia

Osteoporosis

Diagnosis:

Signature _____

Device Details

System S/N.....20052020 Probe S/N.....CMB288
 S/W.....2.1.2260
 RDB.....Caucasian, Female

Sunlight Medical Inc.

Somerset, NJ, 08873, USA

Tel: 732-560-8770 Fax: 732-560-8770

E-mail: Joanne@sunlightnet.com