

# Measurement Report

12/10/01

## Patient Information

First Name.....**Amanda**  
Family Name.....**Johnson**  
Date of Birth.....12/1/01  
Birth Length.....44.7cm  
Birth Weight.....2160.00  
Gestational Age.....32w6d

ID.....2547813  
Mother's First Name.....Rose  
Gender.....Female  
Age.....3w5d



## Measurement Results

Site.....MIDSHAFT TIBIA - CS

SOS.....2920[m/sec]  
Speed of Sound

Z-Score.....-0.7  
Number of standard deviations from the mean of age matched population

Percentile.....24  
Percentile rank to which patient belongs relative to age matched peers



## Findings and Comments

Signature \_\_\_\_\_

## Device Details

System S/N.....20092194    Probe S/N.....CSB105  
S/W.....2.3.2645    Print Date.....12/27/01  
RDB.....Caucasian, Female

## Sunlight Medical Inc

Somerset NJ, 08873, USA  
Tel: 732-560-8770 Fax: 732-560-8771  
E-mail:Joanne@sunlightnet.com